

# Supplementary protocol

---

## NucleoSpin® DNA Stool – Isolation of genomic DNA from chicken feces (Rev. 01)

This supplementary protocol is developed for the isolation of genomic DNA from 200 mg feces from chicken.

*This protocol is only a supplement to the kit's general user manual. Please refer to the kit manual for more detailed information regarding safety instructions, product-specific disclaimers, and especially preparations needed before starting the procedure. The latest version of the user manual is available at [www.mn-net.com/usermanuals](http://www.mn-net.com/usermanuals) or can be requested from our technical service ([tech-bio@mn-net.com](mailto:tech-bio@mn-net.com)). Material safety data sheets (MSDS) can be downloaded from [www.mn-net.com/MSDS](http://www.mn-net.com/MSDS).*

### Additional equipment needed:

- Proteinase K (REF 740506/.30/.75)

### Additional preparations before starting:

- Add Buffer PB to the lyophilized Proteinase K as described in the instruction leaflet.

---

### 1 Prepare sample

Transfer approx. **200 mg chicken feces** to the **NucleoSpin® Bead Tube Type A**.

Add **1 mL Buffer ST1**.

---

### 2 Lyse sample

Vortex **10 min** at **RT** using the MN Bead Tube Holder on Vortex-Genie® 2 at **max. speed**.

Centrifuge for **5 s** to spin down foam.

Add **20 µL Proteinase K** and mix by **shaking** horizontally for **2–3 seconds**.

Incubate at **70 °C** for **30 min**, invert the tube every 10 min to mix the solution.

---

Proceed with step 3 of the standard protocol (Precipitate contaminants).